

Disposable Face Mask (non-Sterile)

Product name: Disposable Face Mask

Size: 17.5 cm x 9.5 cm

Standard: EN 14683:2019 from test report TPHN20031622671

Package quantity: 20 pcs per bag

Structure and composition: Mask sheet, nose clip and mask tape made by non-woven fabric and polypropylene melt-blown fabric.

Application: It is worn by various clinical personnel during non-invasive operation. Covering the user's mouth, nose and jaw, and providing a certain physical barrier to prevent direct transmission of pathogenic microorganisms and particulate matter.

Method: 1. Face Mask nose clip facing outwards, wearing a lanyard.

2. Adjust and press the bridge of the nose bridge to make the mask fit snugly with the nose and cheeks

3. Keep the mask firmly against your face and stretch it so that it completely covers your mouth, nose and chin.

4. When discarding, fold the mask from the inside to the outside, cover the fixing sheet, and bind it with a strao before throwing it into an unrecyclable garbage bin to avoid polluting the environment.

Attention: It is forbidden to use if the packaging is damaged.

This product is a one-time product, and it is strictly forbidden to reuse it.

See the package for the production date and batch number.

Storage: Store in a dry, ventilated, hygienic, non-corrosive gas room.

Images:

